



Comfit™

Comfort by Design.

Disposable and Reusable Rubber Mouthpieces



DIAGNOSTICS

The Spirometry Source

Comfit mouthpieces have been designed to be a low cost alternative to existing rubber mouthpieces and to afford an exceptional mouth seal with a high degree of patient comfort.

Comfit™ disposable mouthpiece

Made in the USA!

Rims dip to accommodate the frenula



Center support prevents air flow restriction

Soft, comfortable material

- Most comfortable for patients
- Lowest cost
- Single patient or reusable models
- Leak-proof seal for better tests
- Fits most PFT valves and filters
- Latex-free

The Comfit is also packaged in kits--ask your dealer for more info!



Comfit-R™ reusable mouthpiece
Same design in autoclaveable Kraton® material!

Your dealer:

**Here is a sample of our end-user literature for our new Comfit Rubber Mouthpieces!
Call SDI at 1-800-678-5782 to request multiple copies--just fill in your information in
the space provided.**

The Comfit will be available in kits with our Pulmoguard and Pulmoguard II filters:

Part No.	Description	Quantity	Dealer Price	List Price
29-7710-025	Comfit Disposable Mouthpiece	25/Box	\$ 21.00	\$ 30.00
29-7710-050	Comfit Disposable Mouthpiece	50/Box	40.50	57.50
29-7720-010	Comfit-R Reusable Mouthpiece	10/Box	13.90	19.90
29-7720-030	Comfit-R Reusable Mouthpiece	30/Box	39.60	56.70
Kits				
29-7927-040	Pulmoguard II filter, Comfit	40/Box	\$ 84.00	\$ 120.00
29-7927-080	Pulmoguard II filter, Comfit	80/Box	154.40	220.00
29-7928-040	Pulmoguard II, Comfit, The Klip noseclip	40/Box	92.40	132.00
29-7928-080	Pulmoguard II, Comfit, The Klip noseclip	80/Box	165.60	236.00
29-7957-040	Pulmoguard filter, Comfit	40/Box	84.00	120.00
29-7957-080	Pulmoguard filter, Comfit	80/Box	154.40	220.00
29-7958-040	Pulmoguard, Comfit, The Klip noseclip	40/Box	92.40	132.00
29-7958-080	Pulmoguard, Comfit, The Klip noseclip	80/Box	165.60	236.00

Call for Quantity Discount Pricing!