

# PRACTICE SAFE SPIROMETRY


## Infection Protection For Your Patient, Your Environment ...and You!

SDI Diagnostics is a leading manufacturer of spirometry supplies with over 20 years experience in PFT filters.

### PulmoGUARD MG™

### MedGraphics® preVent™ users-- protection for you: Save money!

When patients blow during a PFT, millions of potentially harmful pathogens are forcefully blown into your space.



#### Introducing the PulmoGuard MG a member of our new PulmoGuard family

- Feel safer and save money
- Unique, pleated filter media
- Clear plastic for 100% visual inspection
- Recyclable
- Convenient packaging


**sdi**  
DIAGNOSTICS

The Spirometry Source

  
Made In  
America

# PulmoGUARD MG™

- Lower cost
- 99.99% Bacterial and viral efficiency
- Meets ATS standards for air resistance - .862 cm H<sub>2</sub>O/L/Sec at 14/L/Sec
- Dead space 35 ml
- Infection protection for patients and health care professionals



## PulmoGuard MG Reorder Info

Part No.	Filter Description	Qty	Packaging
29-3103-050	PulmoGuard MG filter (50 PulmoGuard MG filters)	50	50
<b>Filter Kits</b>			
29-3103K-050	PulmoGuard MG filter (50 PulmoGuard MG filters, 50 Kushion Klips)	50	50
29-3103C-050	PulmoGuard MG filter (50 PulmoGuard MG filters, 50 Comfit mouthpieces)	50	50
29-3103CK-050	PulmoGuard MG filter (50 PulmoGuard MG filters, 50 Kushion Klips, 50 Comfit rubber mouthpieces)	50	50



Made in America



Your Dealer Is:

