

# PRACTICE SAFE SPIROMETRY

## Infection Protection For Your Patient, Your Environment ...and You!

SDI Diagnostics, a leading manufacturer of spirometry supplies with over 20 years experience in PFT filters.

**PulmoGUARD C™**

**Collins™ users--protection for your patients:  
Save money!**

For your protection, SDI has designed filters specifically for Collins instruments as an alternative to the DC2 filter.



**Introducing the PulmoGuard C  
a member of our new PulmoGuard family**

- Feel safer and save money
- Unique, pleated filter media
- Clear plastic for 100% visual inspection
- Recyclable
- Convenient packaging


**sdi**  
DIAGNOSTICS

The Spirometry Source

  
**Made In  
America**

# PulmoGUARD C™

- Lower cost
- 99.99% Bacterial and viral efficiency
- Meets ATS standards for air resistance - .862 cm H<sub>2</sub>O/L/Sec at 14/L/Sec
- Dead space 69 ml
- Infection protection for patients and health care professionals
- Exceeds all ATS/ERS standards for spirometry



## PulmoGuard C Reorder Info

Part No.	Filter Description	Qty	Packaging
29-3104-050	Pulmoguard C filter (50 PulmoGuard C filters)	50	50
29-3104-100	Pulmoguard C filter (100 PulmoGuard C filters)	100	100
<b>Filter Kits</b>			
29-3104K-050	Pulmoguard C filter (50 PulmoGuard C filters, 50 Kushion Klip noseclips)	50	50
29-3104K-100	Pulmoguard C filter (100 PulmoGuard C filters, 100 Kushion Klips)	100	100
29-3104C-050	Pulmoguard C filter (50 PulmoGuard C filters, 50 Comfit rubber mouthpieces)	50	50
29-3104C-100	Pulmoguard C filter (100 PulmoGuard C filters, 100 Comfit rubber mouthpieces)	100	100
29-3104CK-050	Pulmoguard C filter (50 PulmoGuard C, 50 Kushion Klips, 50 Comfit)	50	50
29-3104CK-100	Pulmoguard C filter (100 PulmoGuard C, 100 Kushion Klips, 100 Comfit)	100	100



Made in America



Your Dealer Is:

