


# Safe-t-CHEK™

## One-Way Valve Mouthpiece for PFT Testing

For asymptomatic low-risk patients only needing expiratory testing, try the Safe-t-chek disposable mouthpiece. Its unique eight-point peripheral support creates a complete seal and its smooth design is easy for patients to use. The Safe-t-chek offers very low air resistance for accurate testing.



The valve closes completely if a patient attempts to inhale. This design virtually eliminates cross-infection: the closed valve interrupts all air flow, preventing the passage of airborne infectious particles.



### Safe-t-chek Fitting List

Manufacturer	Spirometer Model	Fitting
Breon	All	N/N
CDX	All (remove end fitting)	7976
Circadian	All	N/N
Collins	Eagle/Survey	7976
Collins	Survey Plus	N/N
Cosmed	All	N/N
Jaeger	All	7976
KoKo	All	Call
Med. Equip. Designs	Remove Black Fitting	N/N
Micro Direct	All	N/N
Multispiro	All—Remove black fitting	N/N
PK Morgan	232	N/N
Puritan Bennett	Without flow sensor	N/N
Riko	All	N/N
S & M	Portascreen	N/N
SDI	All Astra	7983/7986
SDI	All others	N/N
Sensormedics	922, 2130	7976
Spirolink	All	N/N
Spirometrics	All	N/N
Vacumetrics	All	Call
Vitalograph	All	N/N

### Peak Flow Model/Brand(s)

Air Life	
AsmaPLAN+/Vitalograph	
Assess/Healthscan	
Astech	
Asthma Check/Healthscan	
Boeringer Peak Flow	
Ko-Ko Puf/Spirometrics	
Peak Air/Omron	
The Peak/Baxter, Multispiro	
Personal Best/Healthscan	
Spiroflow/CDX (Same as KoKo Puf)	
Tru-Zone"/Monahagn (remove green fitting)	
Wright, Mini-Wright/Armstrong, Clement-Clarke	

### Cardboard Safe-t-chek /

Plastic Safe-t-chek	Adapter
■	N/N
■	N/N
■	N/N
■	7976
■	N/N
■	7976
■	N/N
■	7979
■	7979
■	N/N

Your dealer:

# sdi

DIAGNOSTICS

The Spirometry Source

800-678-5782

www.sdidiagnostics.com